



... with a focus on
micro level stimulation



"Minus™ Negative Pressure Wound Therapy System A "Revolutionary" Device Offering Better Results

- ✓ Creating a Sub-atmospheric Pressure with Vacuum Source
- ✓ Maintaining Wound Temperature and Moisture
- ✓ Wound Infection and Inflammation

Minus™ Negative Pressure Wound Therapy system orchestrate the wound healing process, Cleaning of wound bed, stabilize the wound bed temperature at optimum, washing the wound bed with a dwelling time to reach dead corners of the wound bed.

Mechanism of Action

Negative pressure interacts at distinct levels with the wound environment. It works on both the lesion's macrostructure and microstructure and combines all these effects into one complicated operation that facilitates wound healing in all wound growth phases. At the macro level, NPWT extracts secretions and fluids from the lesion, decreases periwound edema, increases local blood flow and reduces ulcer size and the possibility of external sources of wound infection.

The stimulation that the negative pressure exerts on the cell surfaces and on many of the cellular elements present in the lesion induces positive changes in both the shape and function of the cells at the micro level. NPWT has been shown to promote angiogenesis, fibrogenesis and the macrophages and leukocytes activity.

A number of NPWT studies have been conducted in virtually all clinical conditions related to wounds and wound management, acute and chronic, from heart surgery to diabetic feet and venous leg ulcers to pressure ulcers. For most of the conditions under which it was tested, its activity has been confirmed.

An interesting choice with NPWT is that it can be integrated into an antiseptic application. This occurs when antiseptics have entered the system to regularly flush wound areas (as was the case with osteomyelitis) or when they are absorbed into the material used to fill the lesion, this occurs. An antimicrobial effect to the indirect debridement of NPWT is added, as the newer antiseptics have a better efficacy and safety profile and have been tested successfully under a range of conditions.

Creating a Sub-atmospheric Pressure with Vacuum Source

The film has an opening where a tube is attached. The tube leads to a vacuum pump and a canister where fluids are collected. The vacuum pump can be set to run or start&stop intermittently using its modes. The vacuum pump removes fluid and infection from the wound. The edges of the wound bring together with the help of vacuum. It also helps the wound healing by promoting the growth of new tissue. When necessary, antibiotics and saline can be introduced into the wound bed.



PRODUCT LIST	REF NO		
NPWT PUMP HMI MAJOR (Non-Sterile, Reusable)	MNSHMIMJR	NPWT FOAM STANDARD MEDIUM (Sterile – Single Use)	MNSFM
NPWT PUMP ADAPTOR (Non-Sterile, Reusable)	MNSADPT	NPWT FOAM STANDARD LARGE (Sterile – Single Use)	MNSFL
NPWT PUMP SERUM HANGER (Non-Sterile, Reusable)	MNSHNG	NPWT FOAM STANDARD EXTRA LARGE (Sterile – Single Use)	MNSFXL
NPWT PUMP CARRYING BAG (Non-Sterile, Reusable)	MNSCBG	NPWT FOAM STANDARD ABDOMINAL (Sterile – Single Use)	MNSFA
NPWT CLOSURE SET O2 DELIVERY EXTRA SMALL (Sterile – Single Use)	MNSO2CSDEXS	NPWT FOAM STANDARD HAND (Sterile – Single Use)	MNSFH
NPWT CLOSURE SET O2 DELIVERY SMALL (Sterile – Single Use)	MNSO2CSDES	NPWT FOAM STANDARD FOOT (Sterile – Single Use)	MNSFF
NPWT CLOSURE SET O2 DELIVERY MEDIUM (Sterile – Single Use)	MNSO2CSDM	NPWT FOAM STANDARD BRIDGE (Sterile – Single Use)	MNSFB
NPWT CLOSURE SET O2 DELIVERY LARGE (Sterile – Single Use)	MNSO2CSDEL	NPWT FOAM SLIM EXTRA SMALL (Sterile – Single Use)	MNSSFXS
NPWT CLOSURE SET O2 DELIVERY EXTRA LARGE (Sterile – Single Use)	MNSO2CSDEXL	NPWT FOAM SLIM SMALL (Sterile – Single Use)	MNSSFS
NPWT CLOSURE SET O2 DELIVERY ABDOMINAL (Sterile – Single Use)	MNSO2CSDEA	NPWT FOAM SLIM MEDIUM (Sterile – Single Use)	MNSSFM
NPWT CLOSURE SET HAND (Sterile – Single Use)	MNSO2CSDEH	NPWT FOAM SLIM LARGE (Sterile – Single Use)	MNSSFL
NPWT CLOSURE SET O2 DELIVERY FOOT (Sterile – Single Use)	MNSO2CSDEF	NPWT FOAM SLIM EXTRA LARGE (Sterile – Single Use)	MNSSFXL
NPWT CLOSURE SET O2 DELIVERY BRIDGE (Sterile – Single Use)	MNSO2CSDEB	NPWT FOAM SLIM ABDOMINAL (Sterile – Single Use)	MNSSFA
NPWT CLOSURE SET EXTRA SMALL (Sterile – Single Use)	MNSXSI	NPWT FOAM SLIM HAND (Sterile – Single Use)	MNSSFH
NPWT CLOSURE SET SMALL (Sterile – Single Use)	MNSSI	NPWT FOAM SLIM FOOT (Sterile – Single Use)	MNSSFF
NPWT CLOSURE SET MEDIUM (Sterile – Single Use)	MNSMI	NPWT FOAM SLIM BRIDGE (Sterile – Single Use)	MNSSFB
NPWT CLOSURE SET LARGE (Sterile – Single Use)	MNSLI	NPWT FILM DRAPE STANDARD (STERIL DETAYI YAZMIYOR)	MNSFDS
NPWT CLOSURE SET EXTRA LARGE (Sterile – Single Use)	MNSXLI	NPWT FILM DRAPE LARGE	MNSFDL
NPWT CLOSURE SET ABDOMINAL (Sterile – Single Use)	MNSAI	NPWT IRRIGATION TUBING N/A (Sterile – Single Use)	MNSTI
NPWT CLOSURE SET HAND (Sterile – Single Use)	MNSHI	NPWT VACUUM ACCESS PORT™ STANDARD (Sterile – Single Use)	MNSVAPS
NPWT CLOSURE SET FOOT (Sterile – Single Use)	MNSFI	NPWT VACUUM ACCESS PORT™ IRRIGATION (Sterile – Single Use)	MNSVAPIR
NPWT CLOSURE SET BRIDGE (Sterile – Single Use)	MNSBI	NPWT VACUUM ACCESS PORT™ O2 DELIVERY (Sterile – Single Use)	MNSVAP02DE
NPWT INCISION DRESSING SET LARGE (Sterile – Single Use)	MNSLID	NPWT VACUUM TUBING N/A (Sterile- Single Use)	MNSVT
NPWT INCISION DRESSING SET MEDIUM (Sterile – Single Use)	MNSMID	NPWT O2 DELIVERY TUBING N/A	MNSO2DET
NPWT INCISION DRESSING SET SMALL (Sterile – Single Use)	MNSSID	NPWT PERILSTATIC PUMP TUBING N/A	MNSPPT
NPWT ZIPPER DRESSING SET LARGE (Sterile – Single Use)	MNSSZD	NPWT Y CONNECTOR TUBING N/A	MNSYT
NPWT ZIPPER DRESSING SET MEDIUM (Sterile – Single Use)	MNSMZD	NPWT EXTRA LONG TUBING N/A (Sterile – Single Use)	MNSXTRLT
NPWT ZIPPER DRESSING SET SMALL (Sterile – Single Use)	MNSLZD	NPWT CONNECTOR N/A (Sterile – Single Use)	MNSCON
NPWT MCT DRESSING SET LARGE (Sterile – Single Use)	MNSSMD	NPWT Y CONNECTOR N/A (Sterile – Single Use)	MNSYCON
NPWT MCT DRESSING SET MEDIUM (Sterile – Single Use)	MNSMMD	NPWT MULTI LUMEN CONNECTOR (Sterile – Single Use)	MNSMLCON
NPWT MCT DRESSING SET SMALL (Sterile – Single Use)	MNSLMD	NPWT CONNECTOR STANDARD (Sterile – Single Use)	MNSCONS
NPWT UNIBODY DRESSING SET XSMALL (Sterile – Single Use)	MNSXSUD	NPWT SERUM TUBING (Sterile – Single Use)	MNST
NPWT UNIBODY DRESSING SET SMALL (Sterile – Single Use)	MNSSUD	NPWT COLLECTION SET 250ML (Non-Sterile)	MNS250
NPWT UNIBODY DRESSING SET MEDIUM (Sterile – Single Use)	MNSMUD	NPWT COLLECTION SET 500ML (Non-Sterile)	MNS500
NPWT UNIBODY DRESSING SET LARGE (Sterile – Single Use)	MNSLUD	NPWT COLLECTION SET 600ML (Non-Sterile)	MNS600
NPWT UNIBODY DRESSING SET XLARGE (Sterile – Single Use)	MNSXLUD	NPWT COLLECTION SET 1000ML (Non-Sterile)	MNS1K
NPWT UNIBODY DRESSING SET SACRUM (Sterile – Single Use)	MNSDS	NPWT CANISTER 250ML N/A (Non-Sterile)	MNSCNS25
NPWT UNIBODY DRESSING SET KNEE (Sterile – Single Use)	MNSDK	NPWT CANISTER 500ML N/A (Non-Sterile)	MNSCNS50
NPWT UNIBODY DRESSING SET CORNER (Sterile – Single Use)	MNSDC	NPWT CANISTER 600ML N/A (Non-Sterile)	MNSCNS60
NPWT UNIBODY DRESSING SET FOOT (Sterile – Single Use)	MNSDF	NPWT CANISTER 1000ML N/A (Non-Sterile)	MNSCNS1K
NPWT UNIBODY DRESSING SET BRIDGE (Sterile – Single Use)	MNSDB	NPWT FILTER CARTRIDGE FOR HMI GAMECHANGER PUMP (Non-Sterile, Reusable)	MNSFLTRGMC
NPWT FOAM STANDARD EXTRA SMALL (Sterile – Single Use)	MNSFXS	NPWT FILTER CARTRIDGE FOR HMI MAJOR PUMP (Non-Sterile, Reusable)	MNSFLTRMJR
NPWT FOAM STANDARD SMALL (Sterile – Single Use)	MNSFS	NPWT FILTER CARTRIDGE FOR 571 PUMP (NonSterile, Reusable)	MNSFLTR571
		NPWT FILTER CARTRIDGE FOR 646 PUMP (NonSterile, Reusable)	MNSFLTR646

Indications

MINUS™ Negative Pressure Wound Therapy System is a non-pharmacological method for the patients who would benefit from a suction device. It is intended to use to physically promote wound healing by means of removal of fluids, including irrigation, and body fluids wound exudates and infectious materials. The patients with the following wounds are indicated: Acute Wounds, Chronic Wounds Pressure Ulcers, Diabetic Foot Ulcers, Venous Leg ulcers, Traumatic wounds, Dehisced wounds, Partial thickness burns, Flaps, and grafts. No Studies are done on the use of Minus™ NPWT System on pregnant and pediatric patients. To use the device on such patients is at surgeons' own disposal. Patient size and age should be considered when prescribing the device. The Minus™ Negative Pressure Wound Therapy System is intended for use in acute and chronic wounds at continual, intermittent and variable mode of settings.

Major HMI pump can be offered with with/without irrigation options.

References

- 1 Lorie, S., Domp Martin, A., Penven, K. et al. Is vacuum-assisted closure a valid technique for debriding chronic leg ulcers? *J Wound Care*. 2004; 13: 249–252.
- 2 Saxena, V., Hwang, C.W., Huang, S. et al. Vacuum-assisted closure: microdeformations of wounds and cell proliferation. *Plast Reconstr Surg*. 2004; 114: 1086–1096.
- 3 Armstrong, D.G., Lavery, L.A., Diabetic Foot Study Consortium. Negative pressure wound therapy after partial diabetic foot amputation: a multicentre, randomised controlled trial. *Lancet*. 2005; 366: 9498, 1704–1710.
- 4 Blume, P.A., Walters, J., Payne, W. et al. Comparison of negative pressure wound therapy using vacuum-assisted closure with advanced moist wound therapy in the treatment of diabetic foot ulcers: a multicenter randomized controlled trial. *Diabetes Care*. 2008; 31: 631–636
- 5 Mouës, C.M., van den Bermd, G.J., Heule, F., Hovius, S.E. Comparing conventional

gauze therapy to vacuum-assisted closure wound therapy: a prospective randomised trial. *J Plast Reconstr Aesthet Surg*. 2007; 60: 672–681.

6 Wanner, M.B., Schwarzl, F., Strub, B. et al. Vacuum-assisted wound closure for cheaper and more comfortable healing of pressure sores: a prospective study. *Scand J Plast Reconstr Surg Hand Surg*. 2003; 37: 28–33.

7 Timmers, M.S., Graafland, N., Bernards, A.T. et al. Negative pressure wound treatment with polyvinyl alcohol foam and polyhexanide antiseptic solution instillation in posttraumatic osteomyelitis. *Wound Repair Regen*. 2009; 17: 278–286.

8 McGuiness, W., Vella, E., & Harrison, D. (2004). Influence of dressing changes on wound temperature. *Journal of Wound Care*, 13(9), 383-385. doi:10.12968/jowc.2004.13.9.26702

9 Wounds UK Best Practice Statement. The use of topical antimicrobial agents in wound management (3rd edn). London: Wounds UK, 2013. Accessed 08.10.2014 at: <http://www.wounds-uk.com/best-practice-statements/best-practice-statement-theuse-of-topical-antimicrobial-agents-inwound-management>



Disclaimer

A surgeon should rely exclusively on his or her own professional medical/clinical judgment when deciding which particular product to use when treating a patient. Maveria Medical Devices Inc. does not dispense medical advice and advocates that surgeons be trained in the use of any particular product before using it in surgery. The content presented is intended to identify the breadth of Maveria Medical Devices Inc. product options. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Maveria Medical Devices Inc. product. Products may not always be available in all markets in a similar capacity due to the local regulatory and/or medical practices. It is advised that you contact a Maveria Medical Devices Inc. representative if you have questions about Maveria Medical Devices Inc. products. Minus™ is a Registered Trademark of Maveria Medical Devices Inc.

Maveria Medical Devices Inc. guarantees that the products are designed, manufactured and packaged on application of the EN ISO 13485 certified quality management with all possible care, using the most suitable state-of-the-art processes and applying the principle of integrating safety into its design and manufacture. This system is designed to guarantee the safe use of device when it is employed under condition and for the purposes for which it is intended, in accordance with precautions described in the operator's manual and/or instructions for use and to minimize the risks associated with the use of the device as far as possible, although they cannot be totally eliminated. The device must only be used under the responsibility of specialist medical personnel, taking account of the unavoidable risks and the possible side effects and complications of the treatment for which, it is designed, including those referred to in other sections of these instructions for use.

U.S. and foreign patents pending

Issue Date Dec 2020

Copyright© 2020. Maveria Medical Devices Inc. All rights reserved.

MNSMJR0620-EN-REV01
Printed in Turkey



MAVERIA MEDICAL DEVICES INC.



Yeşilbayır Mahallesi Sağduyu Caddesi
107/C Mamak/ Ankara-TURKEY

+90 312 341 50 27
+90 312 341 50 26



www.maveramedical.com.tr